

THORAX	ABDOMEN	BONE			
300mA	1/120sec	300mA	1/60sec	300mA	1/15sec
cm	kVp	cm	kVp	cm	kVp
1	60	1	53	1	43
2	62	2	55	2	45
3	64	3	57	3	47
4	66	4	59	4	49
5	68	5	61	5	51
6	70	6	63	6	53
7	72	7	65	7	55
8	74	8	67	8	57
9	76	9	69	9	59
10	78	10	71	10	61
300mA	1/60sec	300mA	1/15sec	300mA	1/4sec
cm	kVp	cm	kVp	cm	kVp
11	80	11	73	11	63
12	82	12	75	12	65
13	84	13	77	13	67
14	86	14	79	14	69
15	89	15	81	15	71
16	92	16	83	16	73
17	95	17	85	17	75
18	98	18	88	18	77
19	101	19	91	19	79
20	105	20	94	20	81
21	109	21	97	21	83
22	113	22	100	22	85
23	117	23	104	23	88
24	120	24	108	24	91
25	120	25	112	25	94
26	120	26	116	26	97
27	120	27	120	27	100
28	120	28	120	28	104
29	120	29	120	29	108